

2.4 Instruments Ordering Information

P/N	Description	Quantity
10160200E	Depth Gauge, Ø2.4	1
10160600E	Drill Guide, Ø0.9/Ø1.7	1
10160700E	Stardrive Screwdriver, Ø1.75	1
10160800E	Screw Holder, Ø2.4	1
10160900E	Drill Bit, Ø1.7	2
10161200E	Canal Cleaning Wire, Ø0.9	1
10160100E	Ø0.9 Guide Wire	5
10160300E	Ø2.5 Cannulated Drill	2
10160400E	Ø1.7 Cannulated Drill Bit	2
10160500E	Ø1.75 Cannulated Stardriver	2
14301700E	Quick Coupling Handle	1
10700200E	Periosteal Elevator, Small	1
10700300E	Retractor, Small	1
10901600E	Screw Forceps	1
10166000E	Screw Box	1
10165000E	Instrument Box	1

3.0 Instruments Ordering Information

P/N	Description	Quantity
10170100E	Screw Holder, Ø3.0	1
10170200E	Ø3.0 Cannulated Drill	2
10170300E	T8 Cannulated Stardriver Shaft	2
14900100E	Ø1.1 Threaded Guide Wire	5
14900200E	Ø1.1 Guide Wire	5
14901201E	Ø2.0 Cannulated Drill with Quick Coupling	2
14900301E	Ø2.0 Cannulated Drill Bit	2
14901400E	Ø2.0 Drill Bit with Quick Coupling, Solid	1
14901500E	Ø2.0 Drill Bit, Solid	1
14301700E	Quick Coupling Handle	1
14900500E	Ø2.0/1.1 Drill Guide	1
14900601E	T8 Stardriver Shaft, Solid	1
14900800E	Ø1.1 Canal Cleaning Wire	1
14900900E	Depth Gauge	1
10700200E	Periosteal Elevator, Small	1
10700300E	Retractor, Small	1
10901600E	Screw Forceps	1
10175000E	Instrument Box	1
10176000E	Screw Box	1

4.0 Instruments Ordering Information

P/N	Description	Quantity
10180400E	Screw Holder, Ø4.0	1
10180600E	Depth Gauge, Ø4.0	1
10180300E	Ø4.2 Cannulated Drill	2
10180500E	Ø2.8 Cannulated Stardriver	2
15000100E	Ø1.25 Threaded Guide Wire	5
15000200E	Ø1.25 Guide Wire	5
15000301E	Ø2.7 Drill Bit with Quick Coupling, Cannulated	2
15000401E	Ø2.7 Cannulated Drill Bit	2
15000601E	Ø4.0 Screw Tap, Cannulated	1
15000900E	Ø2.7/1.25 Drill Guide	1
15000800E	Ø1.25 Canal Cleaning Wire	1
14301700E	Quick Coupling Handle	1
15001401E	2.8 Stardriver Shaft, Solid	1
15001600E	Ø2.7 Drill Bit, Solid	1
15001300E	Ø2.7 Drill Bit with Quick Coupling, Solid	1
10900700E	Screw Forceps	1
10180100E	Instrument Box	1
10180200E	Screw Box	1

4.5 Instruments Ordering Information

P/N	Description	Quantity
10190300E	Canal Cleaning Wire, Ø1.6	1
10190700E	Wire Guide, Ø1.6	1
10190800E	Drill Guide, Ø3.0	1
10191000E	Depth Gauge, Ø4.5	1
10191100E	Screw Holder, Ø4.5	1
10191200E	Tap, Ø4.5, Cannulated	1
10191300E	Drill Bit, Ø3.0	1
10190400E	Guide Wire, Ø1.6	5
10190500E	Drill Bit, Ø3.0, Cannulated	2
10190600E	Drill Bit, Ø4.8, Cannulated	2
10190900E	Ø2.8 Cannulated Stardriver	2
15001401E	Ø2.8 Solid Stardriver	1
10901600E	Screw Forceps	1
14301700E	Quick Coupling Handle	1
10190100E	Instrument Box	1
10190200E	Screw Box	1

5.5 Instruments Ordering Information

P/N	Description	Quantity
10210600E	Wire Guide, Ø1.8	1
10210700E	Drill Guide, Ø3.8	1
10210800E	Depth Gauge, Ø5.5	1
10210900E	Screw Holder, Ø5.5	1
10211000E	Tap, Ø5.5, Cannulated	1
10211100E	Drill Bit, Ø3.8	1
10211200E	Canal Cleaning Wire, Ø1.8	1
10210300E	Ø1.8 Guide Wire	5
10210400E	Drill Bit, Ø3.8, Cannulated	2
10210500E	Drill, Ø5.8, Cannulated	2
10221000E	Ø4.5 Cannulated Stardriver	1
10901600E	Screw Forceps	1
14301700E	Quick Coupling Handle	1
10210100E	Instrument Box	1
10210200E	Screw Box	1

6.5 Instruments Ordering Information

P/N	Description	Quantity
10220600E	Depth Gauge, Ø6.5	1
10220700E	Screw Holder, Ø6.5	1
10220800E	Wire Guide, Ø2.0	1
10220900E	Drill Guide, Ø4.0	1
10221200E	Drill Bit, Ø4.0	1
10221300E	Canal Cleaning Wire, Ø2.0	1
10220300E	Ø2.0 Guide Wire	5
10220400E	Ø4.0 Cannulated Drill Bit	2
10220500E	Drill, Ø6.8, Cannulated	2
10221000E	Ø4.5 Cannulated Stardriver	1
10221100E	Tap, Ø6.5, Cannulated	1
10901600E	Screw Forceps	1
14301700E	Quick Coupling Handle	1
10220100E	Instrument Box	1
10220200E	Screw Box	1

HB Headless Cannulated Compression Screw(II) System

Brochure



TRAUSON®

Trauson(China) Medical Instrument Co., Ltd
 Add:No.9 Long Men Road, Wujin Hi-Tech Industrail Zone, 213164, Changzhou City, Jiangsu Province. P.R.C.
 Tel:86-519-86387075 Fax:86-519-86387279 Email:info@trauson.com http://www.trauson.com

08/13/2006-17.pdf

Products you can trust at a price you can afford

HB Headless Cannulated Compression Screw(II) System

Indications

The HB Headless Cannulated Screw is indicated for: various bone fragment fixation, non-unions and mal-unions.

Please refer to the operative technique for a complete list of indications and contra-indications.

Features & Benefits

- StarDrive recess design provides improved torque transmission.
- Different thread pitch design allows compression in fragments during screw insertion
- Double thread design for easy screw insertion under cortical bone
- Headless design decreases the irritation of tissues and articular surfaces, particularly suitable for periarticular fractures
- Self-tapping groove design
- Deep thread to increase the grasping force
- Two different kinds of thread length - suitable for all types of fractures
- Self tapping drill for convenient, time saving surgery
- Cannulated design allows insertion over guide for exact screw positioning



Implants

HB2.4 Headless Cannulated Screw (II)

P/N(TA)	Screw Length	Threaded Length (L1)
44341010E	10mm	4mm
44341011E	11mm	4mm
44341012E	12mm	4mm
44341013E	13mm	4mm
44341014E	14mm	4mm
44341015E	15mm	4mm
44341016E	16mm	4mm
44341017E	17mm	4mm
44341018E	18mm	4mm
44341019E	19mm	4mm
44341020E	20mm	4mm
44341021E	21mm	4mm
44341022E	22mm	4mm
44341023E	23mm	4mm
44341024E	24mm	5mm
44341025E	25mm	5mm
44341026E	26mm	5mm
44341027E	27mm	6mm
44341028E	28mm	6mm
44341029E	29mm	6mm
44341030E	30mm	7mm
44341032E	32mm	7mm
44341034E	34mm	8mm
44341036E	36mm	9mm
44341038E	38mm	9mm
44341040E	40mm	10mm

HB2.4 Headless Cannulated Screw (II)

P/N(TA)	Screw Length	Threaded Length(L2)
44342016E	16mm	5mm
44342017E	17mm	6mm
44342018E	18mm	6mm
44342019E	19mm	7mm
44342020E	20mm	7mm
44342021E	21mm	8mm
44342022E	22mm	8mm
44342023E	23mm	8mm
44342024E	24mm	8mm
44342025E	25mm	8mm
44342026E	26mm	10mm
44342027E	27mm	10mm
44342028E	28mm	10mm
44342029E	29mm	10mm
44342030E	30mm	12mm
44342032E	32mm	12mm
44342034E	34mm	14mm
44342036E	36mm	14mm
44342038E	38mm	16mm
44342040E	40mm	16mm

HB3.0 Headless Cannulated Screw (II)

P/N(TA)	Screw Length	Threaded Length(L1)
44343010E	10mm	4mm
44343011E	11mm	4mm
44343012E	12mm	4mm
44343013E	13mm	4mm
44343014E	14mm	4mm
44343015E	15mm	4mm
44343016E	16mm	4mm
44343017E	17mm	4mm
44343018E	18mm	4mm
44343019E	19mm	4mm
44343020E	20mm	4mm
44343021E	21mm	4mm
44343022E	22mm	4mm
44343023E	23mm	4mm
44343024E	24mm	5mm
44343025E	25mm	5mm
44343026E	26mm	5mm
44343027E	27mm	6mm
44343028E	28mm	6mm
44343029E	29mm	6mm
44343030E	30mm	7mm
44343032E	32mm	7mm
44343034E	34mm	8mm
44343036E	36mm	9mm
44343038E	38mm	9mm
44343040E	40mm	10mm

HB4.0 Headless Cannulated Screw (II)

P/N(TA)	Screw Length	Threaded Length(L2)
44346020E	20mm	10mm
44346022E	22mm	11mm
44346024E	24mm	12mm
44346026E	26mm	13mm
44346028E	28mm	14mm
44346030E	30mm	15mm
44346032E	32mm	16mm
44346034E	34mm	17mm
44346036E	36mm	18mm
44346038E	38mm	19mm
44346040E	40mm	20mm
44346042E	42mm	21mm
44346044E	44mm	22mm
44346046E	46mm	23mm
44346048E	48mm	24mm
44346050E	50mm	25mm
44346052E	52mm	26mm
44346054E	54mm	27mm
44346056E	56mm	28mm
44346058E	58mm	29mm
44346060E	60mm	30mm
44346064E	64mm	32mm
44346068E	68mm	34mm
44346072E	72mm	36mm

HB3.0 Headless Cannulated Screw (II)

P/N(TA)	Screw Length	Threaded Length(L2)
44344016E	16mm	5mm
44344017E	17mm	6mm
44344018E	18mm	6mm
44344019E	19mm	7mm
44344020E	20mm	7mm
44344021E	21mm	8mm
44344022E	22mm	8mm
44344023E	23mm	8mm
44344024E	24mm	8mm
44344025E	25mm	8mm
44344026E	26mm	10mm
44344027E	27mm	10mm
44344028E	28mm	10mm
44344029E	29mm	10mm
44344030E	30mm	12mm
44344032E	32mm	12mm
44344034E	34mm	14mm
44344036E	36mm	14mm
44344038E	38mm	16mm
44344040E	40mm	16mm

HB4.0 Headless Cannulated Screw (II)

P/N(TA)	Screw Length	Threaded Length(L1)
44345016E	16mm	5mm
44345018E	18mm	6mm
44345020E	20mm	6mm
44345022E	22mm	8mm
44345024E	24mm	8mm
44345026E	26mm	9mm
44345028E	28mm	9mm
44345030E	30mm	10mm
44345032E	32mm	10mm
44345034E	34mm	12mm
44345036E	36mm	12mm
44345038E	38mm	13mm
44345040E	40mm	13mm
44345042E	42mm	14mm
44345044E	44mm	14mm
44345046E	46mm	16mm
44345048E	48mm	16mm
44345050E	50mm	17mm
44345052E	52mm	17mm
44345054E	54mm	18mm
44345056E	56mm	18mm
44345058E	58mm	20mm
44345060E	60mm	20mm
44345062E	64mm	22mm
44345064E	68mm	23mm
44345066E	72mm	24mm

HB4.5 Headless Cannulated Screw (II)

P/N(TA)	Screw Length	Threaded Length
44347025E	25mm	8mm
44347030E	30mm	8mm
44347035E	35mm	8mm
44347040E	40mm	12mm
44347045E	45mm	12mm
44347050E	50mm	12mm
44347055E	55mm	12mm
44347060E	60mm	12mm
44347065E	65mm	12mm
44347070E	70mm	12mm
44347075E	75mm	16mm
44347080E	80mm	16mm
44347085E	85mm	16mm
44347090E	90mm	16mm
44347095E	95mm	16mm
44347100E	100mm	16mm

HB6.5 Headless Cannulated Screw (II)

P/N(TA)	Screw Length	Threaded Length
44349025E	25mm	8mm
44349030E	30mm	8mm
44349035E	35mm	8mm
44349040E	40mm	12mm
44349045E	45mm	12mm
44349050E	50mm	12mm
44349055E	55mm	12mm
44349060E	60mm	12mm
44349065E	65mm	12mm
44349070E	70mm	12mm
44349075E	75mm	16mm
44349080E	80mm	16mm
44349085E	85mm	16mm
44349090E	90mm	16mm
44349095E	95mm	16mm
44349100E	100mm	16mm

HB5.5 Headless Cannulated Screw (II)

P/N(TA)	Screw Length	Threaded Length
44348025E	25mm	8mm
44348030E	30mm	8mm
44348035E	35mm	8mm
44348040E	40mm	12mm
44348045E	45mm	12mm
44348050E	50mm	12mm
44348055E	55mm	12mm
44348060E	60mm	12mm
44348065E	65mm	12mm
44348070E	70mm	12mm
44348075E	75mm	16mm
44348080E	80mm	16mm
44348085E	85mm	16mm
44348090E	90mm	16mm
44348095E	95mm	16mm
44348100E	100mm	16mm